



IPS Implants®

Cranium





Oral and maxillofacial surgery is our passion! We also want to continue our development along with our customers. Day in, day out, we work to develop innovative products and services that satisfy the highest quality demands and contribute to the patient's well-being.

IPS® – Individual Patient Solutions

IPS Implants® Cranium

Traumatic brain injuries, loss of bone integrity or reconstructions as a result of tumors, ulcers or cysts require restoration of the original shape and functions to restore the patient's quality of life.

Despite advances in reconstruction techniques, it has always been a challenge to reconstruct the skull in its original shape.

The use of modern technologies opens up new options in the treatment of complex defect situations. With the development of preoperative virtual planning as well as patient-specific implants, another possibility to achieve a true-to-origin contour of resected bone has been established. IPS® offers matched solutions for the computer-based planning of surgical procedures, the efficient design of customized treatment concepts and the realization of these concepts in the operating theater with functionalized implants and planning aids.



Table of Contents

| | Pages |
|------------------------------------|-------|
| Feature, function and benefit | 6-7 |
| Indications and surgical technique | 8-11 |
| Case studies | 12-13 |
| Osteosynthesis accessories | 14 |
| The IPS® product range | 15 |

Feature, Function and Benefit



IPS® is ideal for solutions customized to the patient by a simple and efficient process – from planning to the functional implant.

With IPS Gate®, we provide a platform which guides surgeons and users reliably and efficiently through the process of inquiring about, planning, and completing patientspecific products. The intuitive concept offers the user maximum mobility, flexibility and functionality. With the "HTTPS" standard, IPS Gate® ensures encrypted data transmission, which is additionally certified by the TÜV Süd seal.

Patient-specific implants, planning aids, and anatomical models are made from various materials using state-of-the-art fabrication technologies. Thanks to computer-based planning and functionalized patient-specific implants, preoperative planning can be implemented in surgery with unprecedented precision.




The resulting advantages for patients are reduced complication rates, improved esthetic and functional results, shortened surgical time and faster rehabilitation.



For more information on the materials
used for IPS Implants® please visit our website

<https://www.klsmartin.com/en/products/individual-patient-solutions-cmf/ips-implants>

IPS Implants® Cranium

| | Feature and function | Benefit |
|--------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Planning process  | <ul style="list-style-type: none"> Simple and efficient interaction with the user via IPS Gate® Planning, fabrication, shipping from a single source Provision of three-dimensional planning data Planning time 5-7 working days | <ul style="list-style-type: none"> Maximum mobility, flexibility and functionality Complete service with the requirement for coordinating multiple services eliminated High degree of safety in planning Save time with efficient case processing |
| Marking guides  | <ul style="list-style-type: none"> Enable transfer of virtual planning to the OR In cranial reconstruction: enable precise determination of the implant position For the correction of congenital malformations: integration of bone rearrangement in one or more marking guides Manufactured from polyamide | <ul style="list-style-type: none"> High planning and implementation reliability Exact fit of the implant on the defective area Positioning aid for the correct arrangement of the bone segments High biocompatibility |
| Implants  | <ul style="list-style-type: none"> High variety of materials Option of overlapping the defect or covering it with a perfect fit Latest production technologies such as additive manufacturing marPOR (UHMWPE, ultra-high-molecular-weight polyethylene) <ul style="list-style-type: none"> Porous 3D interconnecting pore structure Delivered sterile High-performance polymer PEEK (polyether ether ketone) Titanium mesh and solid titanium Implant based on the individual patient CT scan, checked for perfect fit ex-works | <ul style="list-style-type: none"> Wide range of choices within the context of best possible patient care Complete freedom of design for implants Maximum flexibility and stability Allowing soft tissue ingrowth <ul style="list-style-type: none"> Flexible and strong material Good drainage properties Can promote cell ingrowth, vascularization and osteointegration Ready-to-use Physical properties similar to those of human cortical bone Intraoperative adjustment possible if required High implant stability Cost-effective alternative through patient-specific preforming of the implant No sharp edges, as cutting or bending to size is no longer necessary |

Step by Step to Optimal Treatment

Indications

Cranial reconstructions due to infections, tumors, cysts or rejection reactions.

Restoration of form and function following trauma.

Correction of congenital deformities (craniofacial malformations).



- Cranial reconstruction with
- Additively manufactured titanium implant (AMTi)
 - marPOR implant (UHMWPE)
 - PEEK implant
 - Standard titanium mesh
 - Standard solid titanium

(optional: use of a marking guide)



Correction of craniofacial malformations using marking guides



Surgical Technique

Cranial reconstruction with PEEK implant and marking guide

Pages 10-11





Virtual planning

To create the case the patient data and other case-related information are uploaded to the web-based platform IPS Gate®.

The data is prepared for case planning on the basis of the user's requirements and information. An integrated chat function and web meetings are available for direct communication between the IPS® developer and user.

In consultation with the user, the defective regions are identified and the resection limits defined.

Then the marking guide and a case-specific optimized implant are created. The type, diameter and size of the implant are defined exactly according to the user's wishes and are based on anatomical and clinical conditions.

At the end the user approves the design for production.

Note:

More detailed information on how to prepare patient data for virtual planning can be found in our brochure "IPS Implants® Scan Protocol Cranium / Midface / Midface Orbita / Mandible / Mandible Reconstruction".


Resection of the defective region

To prepare the cranial vault, the marking guide is placed and fixed.

This specifies the resection lines and angles for the preparation and resection of the defect. For example, marking can be performed with a piezo device.

Then the marking guide is removed and the resection performed along the marked line.





Placement of the implant

After the defective region has been prepared using a marking guide, the IPS® implant is placed as a next step.

Based on the information stored in the marking guide, a fit with millimeter precision can be guaranteed.



Fixation of the implant

The IPS® implant is fixated to the native bone with osteosynthesis screws (e.g. 1.5-mm maxDrive® Drill-Free screws).

PEEK implants can also be fixated with osteosynthesis plates (e.g. with 1.5-mm Low Profile Neuro System).

Note:

In addition to the IPS® implant and the included marking guides, the required osteosynthesis accessories (KLS Martin osteosynthesis screws and plates in the planned diameters and lengths as well as the corresponding screw-driver and, if applicable, twist drill) must be available in sterile condition. They are not included in the IPS® package.



Cranial reconstruction with additively manufactured titanium implant
Restoration with laser-sintered titanium implant with osteo-conductive mesh structure



Cranial reconstruction with titanium mesh
Restoration with patient-specific preformed titanium mesh



Cranial reconstruction with solid titanium
Restoration with patient-specific preformed solid titanium



Cranial reconstruction with PEEK implant
Restoration with perforated PEEK implant to enable connective tissue ingrowth



Cranial reconstruction with marPOR implant (UHMWPE)
Restoration with porous polyethylene implant to enable osteointegration, vascularization and connective tissue ingrowth



Cranial reconstruction with a solid PEEK implant
Restoration with unperforated PEEK implant with marking guide for precise resection of the defect



Correction of craniofacial malformations
Restoration using additively manufactured marking guides made of polyamide


Osteosynthesis Accessories



In addition to the IPS® implant and/or the included marking guides, the following osteosynthesis accessories are required:

- Sufficient number of KLS Martin osteosynthesis screws in the planned diameters (1.5 mm) for fixating the implant and marking guides.
- When using a marPOR (UHMWPE) or PEEK implant: sufficient number of KLS Martin osteosynthesis plates (L1® Cranium or L1® Midface).
- A screwdriver to fit the planned osteosynthesis screws.
- If no Drill-Free screws are used: a twist drill suitable for the planned osteosynthesis screws.

The IPS® Product Range



IPS CaseDesigner®

The IPS CaseDesigner® makes virtual 3D surgical planning easier and faster than ever before. With this flexible software tool, orthognathic procedures can be efficiently and reliably planned and simulated, and then applied to treatment in the operation in a customized manner.



IPS Gate®

The web-based platform and app guide surgeons and users reliably and efficiently through the process of inquiring about, planning, and completing patient-specific products. With the HTTPS standard IPS Gate® guarantees encrypted data transmission, which is additionally certified by the TÜV Süd seal.



IPS Implants®

Patient-specific implants, planning aids, and anatomical models are made from various materials using state-of-the-art fabrication technologies. Thanks to computer-based planning and functionalized patient-specific implants, preoperative planning can be implemented in surgery with unprecedented precision.



KLS Martin Group

KLS Martin Australia Pty Ltd.

Sydney · Australia
Tel. +61 2 9439 5316
australia@klsmartin.com

KLS Martin do Brasil Ltda.

São Paulo · Brazil
Tel. +55 11 3554 2299
brazil@klsmartin.com

KLS Martin Medical (Shanghai) International Trading Co., Ltd

Shanghai · China
Tel. +86 21 5820 6251
info@klsmartin.com

KLS Martin India Pvt Ltd.

Chennai · India
Tel. +91 44 66 442 300
india@klsmartin.com

KLS Martin Italia S.r.l.

Milan · Italy
Tel. +39 039 605 67 31
info@klsmartin.com

KLS Martin Japan K.K.

Tokyo · Japan
Tel. +81 3 3814 1431
info@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Penang · Malaysia
Tel. +604 261 7060
malaysia@klsmartin.com

KLS Martin de México, S.A. de C.V.

Mexico City · Mexico
Tel. +52 55 7572 0944
mexico@klsmartin.com

KLS Martin Nederland B.V.

Huizen · Netherlands
Tel. +31 35 523 45 38
infonl@klsmartin.com

KLS Martin SE & Co. KG

Moscow · Russia
Tel. +7 499 792 76 19
russia@klsmartin.com

KLS Martin Taiwan Ltd.

Taipei · Taiwan
Tel. +886 2 2325 3169
taiwan@klsmartin.com

KLS Martin SE & Co. KG

Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

KLS Martin UK Ltd.

Reading · United Kingdom
Tel. +44 118 467 1500
info.uk@klsmartin.com

KLS Martin LP

Jacksonville · Florida, USA
Tel. +1 904 641 77 46
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.

Hanoi · Vietnam
Tel. +49 7461 706-0
info@klsmartin.com

KLS Martin SE & Co. KG

A company of the KLS Martin Group

KLS Martin Platz 1 · 78532 Tuttlingen · Germany
P.O. Box 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com